ALLERGAN

Marlow International, Parkway, Marlow,
Bucks, SL7 1YL, United Kingdom.

T. +44 (0)1628 494456

F. +44 (0)1628 494956

E. Productsupport_eame@allergan.com

www.allergan.com

NATRELLE® and ALLERGAN® and the Representation
of an Eye Logo are registered trademarks of
Allergan, Inc. All other trademarks mentioned herein
are property of their respective owners.

© 2015 Allergan. All rights reserved.

3441-01 (L3441 Rev.03) 08/2015

- EN BREAST IMPLANTS AND TISSUE EXPANDERS
- FR IMPLANTS MAMMAIRES ET EXPANSEURS TISSULAIRES
- DE BRUSTIMPLANTATE UND GEWEBEEXPANDER
- IT PROTESI MAMMARI ED ESPANSORI TISSUTALI
- ES IMPLANTES MAMARIOS Y EXPANSORES TISULARES
- PT IMPLANTES MAMÁRIOS E EXPANSORES TECIDULARES
- NL BORSTIMPLANTATEN EN WEEFSELEXPANDERS
- SV BRÖSTIMPLANTAT OCH VÄVNADSEXPANDERS
- NO BRYSTIMPLANTATER OG VEVSEKSPANDERE
- DA BRYSTIMPLANTATER OG VÆVSEKSPANDERE
- FI RINTAIMPLANTIT JA -KUDOSVENYTTIMET
- EL ΕΜΦΥΤΕΥΜΑΤΑ ΜΑΣΤΟΥ ΚΑΙ ΔΙΑΣΤΟΛΕΙΣ ΙΣΤΩΝ
- TR MEME İMPLANTLARI VE DOKU GENİŞLETİCİLERİ
- RU ГРУДНЫЕ ИМПЛАНТАТЫ И ТКАНЕВЫЕ ЭКСПАНДЕРЫ
- UK ІМПЛАНТАТИ ТА ТКАНИННІ ЕКСПАНДЕРИ МОЛОЧНОЇ ЗАЛОЗИ

NATRELLE® BREAST IMPLANTS AND TISSUE EXPANDERS

DESCRIPTION

Natrelle® breast implants and tissue expanders are designed for use in augmentation and reconstruction mammoplasty. All Natrelle® implants and expanders are constructed of a silicone elastomer shell and are latex free.

- Natrelle® single lumen gel-filled breast implants are pre-filled with cohesive silicone gel designed to simulate natural breast tissue.
- Natrelle® saline-filled breast implants are filled with saline at the time of surgery.
- Natrelle® 133 tissue expanders are intended for temporary subcutaneous implantation to develop surgical flaps and additional tissue coverage.
- Natrelle® 150 double lumen gel/saline breast implants are designed to function as both tissue expanders and long-term breast implants for one-stage breast reconstruction or augmentation.
- Natrelle® breast implants and tissue expanders contain no latex or natural rubber materials.

IMPLANT DESIGN FEATURES

Gel and Gel/Saline Implants

- All Natrelle® textured gel and gel/saline implants have a BIOCELL™ textured surface engineered with a deep open pore design for firm tissue adherence.
- The INTRASHIEL™ shell features a patented barrier coat between two layers of silicone elastomer to minimise gel diffusion.
- The Natrelle® 150 includes a magnetic resonance imaging (MRI) compatible self-sealing Mini Remote Injection Site which contains a titanium needle guard to prevent inadvertent puncture through the base of the injection site.
- Anatomically shaped Natrelle® single lumen gel-filled breast implants include orientation dots to aid in correct implant positioning during surgery (see Figure 1).


FIGURE 1

Saline Implants

- Diaphragm valves in Natrelle® saline-filled breast implants are designed for ease in filling and subsequent air removal.

Tissue Expanders

- BIOCELL™ textured surface is designed to promote tissue adherence.
- The stable base in the Natrelle® 133 provides greater control over expansion direction.
- Natrelle® 133 tissue expanders, with integral MAGNA-SITE™ injection site are supplied with a MAGNA-FINDER™ external locating device. The MAGNA-SITE™ and MAGNA-FINDER™ contain rare-earth, permanent magnets for an accurate injection site locating system. In vitro tests show that the MAGNA-SITE™ is detectable through 60mm of phantom tissue.
- All injection sites are self-sealing and contain a titanium needle guard to prevent inadvertent puncture through the base of the injection site (see Figure 2).


FIGURE 2

MAGNA-SITE™ & MAGNA-FINDER™ Locating System

NATRELLE® ACCESSORIES

- Fill Tube Plug Kits;
- MAGNA-FINDER™
- Other product accessories are available separately.

EACH PATIENT MUST BE INDIVIDUALLY EVALUATED FOR IMPLANT SURGERY BASED ON THE CLINICAL JUDGEMENT OF A QUALIFIED SURGEON.

INDICATIONS

- Unilateral or bilateral hypoplasia of the breast.
- Breast reconstruction in patients with adequate tissue covering following mastectomy or trauma.
- Asymmetry, ptosis, or aplasia of the breast.
- Replacement of implants for medical or cosmetic reasons.
- Congenital deformity of the breast.
- Breast reconstruction in patients following mastectomy or trauma in the case of tissue expanders.
- Treatment of soft tissue deformities in the case of tissue expanders.
- A patient deemed suitable for breast augmentation must be at least 18 years old (22 years old in Singapore).

CONTRAINDICATIONS

- Tissue covering determined inadequate or unsuitable by the surgeon.
- Active infection, local and systemic.
- Existing carcinoma of the breast without mastectomy and residual gross local tumour of the breast after mastectomy.
- Advanced fibrocystic disease considered to be pre-malignant without mastectomy.
- Use of drugs that might result in high surgical risk and/or significant postoperative complications, including drugs that would interfere with blood clotting.
- A patient that demonstrates or shows signs of psychological instability (i.e., an inappropriate attitude or motivation).
- Women who are currently pregnant or breastfeeding.
- Natrelle® 133 tissue expanders contain a MAGNA-SITE™ and should not be used in patients who already have implanted devices that would be affected by a magnetic field (e.g., pacemakers, drug infusion devices).
- Diagnostic testing with MRI is contraindicated in patients with Natrelle® 133 tissue expanders in place. The MRI equipment could cause movement of the MAGNA-SITE™ tissue expander and result in not only patient discomfort but also tissue expander displacement, requiring revision surgery. Also, the MAGNA-SITE™ magnet can interfere with MRI and X-ray detection capabilities. All other Natrelle® implants are MRI safe.

RELATIVE CONTRAINDICATIONS/PRECAUTIONS

- Ptotic breasts where nipple falls below the inframammary fold, without concurrent mastopexy.
- To varying degrees, radiation damage, ulceration, compromised vascularity, or history of compromised wound healing which may affect tissue covering suitability.
- Previous repeated contour correction failures.
- Patients about to undergo radiation therapy and/or chemotherapy as this may make the use of breast implants and tissue expanders more difficult and increase the risk of complications.
- Physiological condition determined by the surgeon to pose unduly high risk of surgical and/or postoperative complications. To varying degrees, obesity, smoking, diabetes, autoimmune disease, coagulopathy, chronic lung or severe cardiovascular disease may affect patient suitability for surgical implantation.

INFORMATION THAT SHOULD BE PROVIDED TO THE PATIENT

All patients should be informed of all the potential benefits and risks (see WARNINGS below) associated with the procedure prior to surgery.

Expected benefits include facilitating emotional healing after cancer, eliminating external prostheses, regaining body symmetry, allowing freedom in clothing and physical activities, and improving sexual or interpersonal relationships.

Patients should be informed about the implant options available, the surgical procedure including implant placement and incision site options. As this surgery will most likely be carried out under a general anaesthetic patients should be made aware of the risks associated with anaesthesia. Patients should discuss with their surgeon any history which may indicate a contraindication (relative or absolute) to surgery. Post-surgery care should be discussed, including time for wound healing, the need for any drainage tubes, recuperation duration and the need for implant integrity to be evaluated on a regular basis after primary surgery. Patients should be advised to consult a physician or pharmacist before using topical medicines (e.g. steroids) in the breast area, and if any clinical examination or surgery in the breast area is planned the patient should inform the doctor or nurse of the presence of an implant. The surgeon should advise the patient to consult a physician should she suspect any complications. All patients should receive a patient information booklet provided by Allergan.

Once the patient has received all the information, she should take at least 30 days to think about the risks and benefits of having breast implants before making a final decision.

WARNINGS

The surgeon should advise the patient that management of the complications listed below may include additional surgery or explantation. Tissue expander patients should be advised that tissue expanders are only to be used for a short term until the tissue has expanded sufficiently.

Breast implants have a limited lifetime and the implant may have to be removed or replaced which may necessitate revision surgery. Various factors, including the type of implant inserted, the type of surgery, injury to the breast, and excessive repetitive compression of the implant, may impact the longevity of the implants. Details on the expected lifetime of the implants are presented in the Rupture/Deflation section below. As many factors affect the lifetime of a device and are outside the control of the manufacturer, the life expectancy of the implant cannot be guaranteed. The surgeon should discuss the necessity of pre-screening mammography with each patient as appropriate for her age and medical history.

1. Rupture/Deflation

Gel implants may rupture, and saline or gel/saline implants may deflate at any time and require replacement or revision surgery. As ruptures are most often clinically silent, a radiological assessment may be required to aid diagnosis.

Causes of rupture or deflation include:

- Damage by surgical instruments; puncture of the valve may also occur from improper insertion of the fill tube for saline implants.
- Other trauma during surgery, such as improper handling or manipulation.
- For *Natrelle*® 150 double lumen implants, *Natrelle*® saline implants and *Natrelle*® 133 tissue expanders, underfilling below the recommended fill volume range for the individual implant may result in folds, fold abrasion, and potentially, crease-fold failure; overfilling above the recommended fill volume range for the individual implant may compromise shell integrity.
- Capsular contracture, or abrasive calcifications in the fibrous capsule.
- Closed or external capsulotomy.
- Stressors such as trauma, intense physical activity, vigorous massage and/or manipulation.
- Excessive compression during mammographic imaging.
- Leakage through remote port devices or through an unsealed or damaged valve.
- Umbilical endoscopic-assisted approach; preliminary reports indicate that there may be a higher incidence of deflation with this approach.

Long term Allergan Post-Market Surveillance data over fourteen years on single lumen and double lumen gel/saline breast implants indicates a rupture rate between 0.37%-1.09%. Allergan US clinical study data on gel implants indicates a rupture rate between 7.7% -9.7% at 10 years.

Long term Allergan Post-Market Surveillance data over fourteen years on single lumen saline-filled breast implants and tissue expanders indicates deflation rates between 0.267%-6.99%. Published results from Allergan US clinical study data indicated a deflation rate of 10.5% for saline implants at 10 years.

2. Capsular Contracture

Formation of a fibrous tissue capsule around an implanted device is a normal physiological response. Fibrous capsular contracture remains a common complication following breast implant surgery and is one of the most common reasons for reoperation. The cause of capsular contracture is unknown, however it is most likely multifactorial and may be more common following infection, haematoma, and seroma. Contracture develops to varying degrees, unilaterally or bilaterally, and may occur within weeks to years after surgery. Contracture of the fibrous capsular tissue surrounding the implant may cause a range of symptoms including firmness, discomfort, pain, distortion, palpability, and/or displacement. Severe cases are considered the most clinically significant, and may require surgical intervention. Capsular contracture may recur subsequent to corrective surgical procedures. DO NOT treat capsular contracture by external compression or massage, which may result in implant damage, deflation, folds, and/or haematoma.

3. Infection

Infection around a breast implant may occur within days, weeks, or even years, after surgery. Signs of acute infection reported in association with implants include erythema, tenderness, fluid accumulation, pain, and fever. Erythema may also occur as a normal response to expansion. Infection that is unresponsive to treatment may require implant removal. Very rarely, Toxic Shock Syndrome has been reported as a possible complication of breast implant surgery and may also be associated with other types of implant surgery.

4. Necrosis

Necrosis may inhibit wound healing and require surgical correction and/or explantation. Permanent scar deformity may occur as a result of necrosis. Placement, expansion and pressure of the remote injection site (in the *Natrelle*® 150) may induce necrosis particularly with unsuitable skin flaps. Do not use microwave diathermy in patients with breast implants. Microwave diathermy has been reported to cause tissue necrosis, skin erosion, and implant extrusion.

5. Haematoma/Seroma

Haematoma/seroma may occur in the postoperative period inhibiting wound healing, or have delayed onset, either of which may require surgical correction and/or explantation.

6. Inflammatory Reaction

Studies evaluating the capsules around textured tissue expanders have reported possible silicone particles within giant cells, indicative of a local (and non-specific) foreign body reaction, and silicone granuloma formation. Another study suggests that certain types of capsule cells, including some perceived as giant cells, may actually be secretory cells that form in response to the frictional forces of the tissue expander, providing lubrication at the capsule-expander interface. In case of an inflammatory reaction, the surgeon is advised to remove the device from the patient's body and to secure any evidence on the possible cause of the inflammatory reaction and treat the patient correspondingly. It is advised not to replace the implant until the inflammatory reaction has passed completely and its cause has been eliminated.

7. Extrusion

Unstable or compromised tissue covering and/or interruption of wound healing may result in extrusion of the implant. In case of an extrusion, the device should be regarded as contaminated and should be removed. It may be replaced with another device after the wound has sufficiently healed.

8. Wrinkling and Folds

Palpable, or even visible, wrinkles and folds may occur (this may be caused by underfilling in the case of saline implants). Folds may result in thinning and erosion of adjacent tissue, and extrusion of the implant. Folds may also result in crease-fold failure and rupture/deflation of the implant. If wrinkling occurs, the device may be replaced with an implant with a different filler or shape.

9. Interference with Standard Mammography/Self-Examination

The patient should continue to perform regular breast examinations for cancer screening; however, this may be more difficult with an implant.

The patient should be informed by the physician about the possible interference of the implant on the self-examination of the breast.

Patients should be instructed to inform their radiologists of the presence of an implant. With breast implants, routine screening mammography will be more difficult as the implant may interfere with diagnostic imaging. Because the breast and implant are squeezed during mammography, an implant may rupture during the procedure. More x-ray views may be necessary for women with breast implants; therefore, a patient may receive more exposure to radiation. However, the benefit of mammography is likely to outweigh the risk of the additional x-rays. Ultrasound may be a useful adjunct to mammography. Breast tissue imaging may be improved by submuscular placement of the implant. Orientation marks on anatomical gel-filled implants may be visible on mammographic images.

10. Pain

As expected following any invasive surgical procedure, pain of varying intensity and duration may occur following implantation. In addition, improper size, placement, surgical technique, or capsular contracture may result in pain associated with nerve entrapment or interference with muscle motion. Unexplained pain must be promptly investigated.

In the case of tissue expanders, the expansion process may cause some discomfort, but should not cause excessive pain. Pain may indicate expansion beyond tissue tolerance, which could result in ischemia and necrosis. Further expansion should be discontinued until the pain is resolved.

11. Breast Feeding and Sensation

Sensation in the nipple and breast can increase or decrease after implant surgery, is typically lost after complete mastectomy where the nipple itself is removed, and can be severely lessened by partial mastectomy. Breast implants may impact the ability to breast feed, though there is no conclusive clinical study data to support this. The periareolar incision may be associated with a

higher likelihood of breast feeding difficulties than other incision sites. The risk of temporary or permanent changes in breast sensation resulting from breast surgery could interfere with the patient's ability to breast feed. Nerve traction and compression have been reported in rare cases in association with tissue expansion. For saline implants, immediate partial deflation should be a standard precaution if nerve impingement is suspected, and filling should not resume until the problem is resolved.

12. Dissatisfaction with Cosmetic Results

Scar deformity, hypertrophic scarring, capsular contracture, asymmetry, displacement, incorrect size, unanticipated contour, palpability, and sloshing (*Natrelle*® 150 and saline-filled implants), may occur. In some cases, cosmetic concerns may also lead to medical concerns. Careful surgical planning and technique can minimise, but not preclude, the risk of such results. Pre-existing asymmetry may not be entirely correctable. Revision surgery may be indicated to maintain patient satisfaction but carries additional considerations and risks. If the patient is dissatisfied with the cosmetic result, revision surgery may be indicated; the device can be replaced with another device of different height, width, projection, volume, shape or filling, or may be placed in a different position in order to achieve a cosmetic result which is more pleasing to the patient.

Re-positioning of the implant during subsequent procedures should be carefully evaluated by the medical team and care taken to avoid contamination of the implant. Use of excessive force during any subsequent procedure can contribute to localized weakening of the breast implant shell potentially leading to decreased device performance.

13. Rotation

Rotation of an anatomical implant may occur. Proper placement and pocket dissection reduces the risk of occurrence. Revision surgery may be necessary to correct rotation. In case of rotation, it is advised to rotate the device back into its correct position in an open surgical procedure. Reshaping of the implant pocket may be necessary to avoid any further rotation in the future.

14. Ptosis

Ptosis occurs naturally in all breasts over time. In case of ptosis, a mastopexy may be performed and/or the device may be replaced by another product with a larger volume or greater projection.

15. Distortion

Tissue expansion is a time and labour intensive process that may cause temporary discomfort and distortion. If distortion occurs, the cause should be analysed and eliminated.

16. Calcification

Calcification commonly occurs in mature breast tissue with or without implantation. Microcalcification after implantation typically occurs on or around the fibrous capsule in thin plaques or accumulations. Extensive microcalcification may cause breast hardness and discomfort, and may necessitate surgical intervention.

17. Tissue Atrophy/Chest Wall Deformity

Pressure of a breast implant or expander may cause tissue atrophy. In rare cases chest wall deformity has also been reported in association with the use of breast implants and tissue expanders.

18. Gel Diffusion

Minute quantities of silicone may diffuse through the elastomer envelope of gel-filled implants. The detection of small quantities of silicone in the periprosthetic capsule, axillary lymph nodes, and other distal regions in patients with apparently unruptured, conventional gel-filled implants has been reported in the medical literature. However, there has been only limited evidence in medical literature associating gel diffusion with local complications in breast implant patients. If significant gel diffusion occurs, the device should be checked for any possible leakage or flaws.

19. Adulterated Fill

Do not use adulterated fill. Saline implants and tissue expanders are to be filled only with sterile saline for injection, and only as described in Instructions for Use.

20. Inadequate Tissue Flap

Inadequate tissue flap following expansion may occur and may require additional surgery and expansion.

21. Deformation

The unique nature of the highly cohesive silicone implant may require a larger incision compared to the incision size required for other silicone-filled implants to avoid skin edge trauma, implant deformation or separation/disruption of the gel. Excessive force upon insertion of the implant may compromise the precisely defined shape of the device, potentially leading to an undesirable cosmetic outcome.

RESEARCH ON SILICONE IMPLANTS

A report published in 1998 by a US National Science Panel, appointed by Judge Sam Pointer, evaluated the scientific data on silicone breast implants in relation to connective tissue diseases and immunologic dysfunction. No association was found between silicone gel-filled implants and any of the definite connective tissue disorders (including Sjögren's Syndrome) or other autoimmune/rheumatic conditions. They found that women with silicone breast implants do not display a silicone-induced systemic abnormality in the types or functions of cells of the immune system.

In 1999, an independent review from a committee at the Institute of Medicine in the US reported that connective tissue disorders, cancer, neurological diseases or other systemic complaints or conditions are no more common in women with breast implants than in women without implants. They concluded that a review of the toxicology studies of silicones and other substances known to be in breast implants does not provide a basis for health concerns.

Lymphoma, including anaplastic large T-cell lymphoma (ALCL) – Information from medical literature has suggested a possible association, without evidence of causation, between breast implants and the very rare occurrence of ALCL in the breast. The disease is exceptionally rare, may present as a late occurring peri-prosthetic seroma, and occurs in women with and without breast implants. Specific testing is needed to distinguish ALCL from breast cancer. The majority of the reported cases had an indolent clinical course following capsulectomy with or without adjuvant therapy, which is generally uncharacteristic of systemic ALCL. Treatment should be determined in consultation with a hemato-oncologist.

INSTRUCTIONS FOR USE**SURGICAL PROCEDURE**

Allergan relies on the surgeon to know and follow the proper surgical procedures with *Natrelle*® implants. The surgeon must carefully evaluate implant size and contour, incision placement, pocket dissection and implant placement criteria with respect to the patient's anatomy and desired physical outcome. Planning should include clear delineation of aesthetic goals to ensure mutual understanding between surgeon and patient. The surgeon should observe current and accepted techniques to minimise the risk of adverse and potentially disfiguring reactions.

Natrelle® products are designed and tested for compatibility with sterile water and saline solution. Other substances, such as alcohol or other chemical agents have not been tested in combination with the *Natrelle*® products.

Do not immerse the implant in povidone-iodine solution (for example, *Betadine*®). If this solution is used in the pocket, ensure that it is rinsed thoroughly so no residual solution remains in the pocket

Natrelle® products should not be exposed to extreme heat, cold or pressure. No excessive force should be used while implanting or removing an implant or expander, and, accordingly, the skin incision should be planned for a sufficient size. No sharp objects such as knives or needles should be used in direct vicinity of any *Natrelle*® product, with the exception of the instruments used to fill an expander or saline-filled implant via the specially designed valve or port.

Do not alter the implants or attempt to repair or insert a damaged device.

Do not place more than one implant per breast pocket.

Back-up implants must be available during the procedure.

SINGLE USE

These products are intended for single use only.

DO NOT reuse explanted products.

DO NOT attempt to re-inflate saline implants following implantation.

RISKS ASSOCIATED WITH REUSE

Natrelle® breast implants and tissue expanders are not intended to be re-sterilized or re-used. The cleaning and autoclaving process can cause damage to the breast implants/tissue expanders, which could lead to loss of structural integrity. Reuse of the device can cause risk of infection to the patient.

PRODUCT IDENTIFICATION

Product labels are supplied within the internal product packaging of each *Natrelle*® implant. The product labels provide specific information which allows product identification.

Important: These labels must be attached to the patient and hospital/doctors records to ensure product identification and device traceability.

STERILE PRODUCT

Each implant is sterilised by dry heat sterilisation and is supplied in a sealed, double primary package.

STORAGE CONDITIONS

Avoid prolonged exposure to extreme storage conditions. Store these devices at ambient room temperatures and at atmospheric pressure and in dry conditions away from direct sunlight.

HOW TO OPEN STERILE PRODUCT PACKAGE

Remove the implant and accessories (where applicable) from their packages in an aseptic environment and using talc-free gloved hands.

DO NOT expose the implant to lint, talc, sponge, towel, skin oils or other surface contaminants.

1. A non-sterile team member peels open the outer package.
2. The surgeon/scrub nurse removes the inner package and places it into the sterile field.
3. Peel open the inner package.
4. Gently retrieve the implant.

Prior to use, keep the implant covered in the inner package, to prevent contact with airborne and surgical field particulate contaminants.

PRELIMINARY PRODUCT EXAMINATION

Prior to use, examine the implant for any evidence of damage or particulate contamination.

GEL AND GEL/SALINE IMPLANTS**IMPLANT PLACEMENT**

Ensure incision is sufficiently large, particularly for textured implants, to facilitate insertion and avoid damage to the device. Inadequate pocket dissection increases the risk of rupture and implant malposition. DO NOT use excessive force during placement of gel-filled implants. Silicone gel may be permanently deformed due to over-manipulation, resulting in deformation of the shape.

NATRELLE® 150 IMPLANT PLACEMENT AND FILLING TECHNIQUES**1. Expandable Breast Implant Placement**

Plan and dissect the surgical pockets for placement of the implant and the Mini Remote Injection Site, using current and accepted surgical techniques. Precise pocket dissection is recommended, and over-dissection should be avoided.

- a) Place the implant flat and correctly oriented in the pocket.
- b) Place the Mini Remote Injection Site flat and correctly-oriented in a separate subcutaneous pocket, ensuring its palpability. The site for the

remote valve pocket should be generally positioned lateral (not inferior) to the Breast in the mid-axillary line. Subcutaneous sutures may be used to maintain the position of the pocket.

- c) Ensure that the injection site connector tubing is not kinked, as this may restrict fluid flow during expansion.

2. Remote Injection Site Preparation

Always verify the location and orientation of the injection site prior to filling. Ensure that the base of the injection site is properly oriented for needle entry.

3. Expandable Mammary Implant Filling

If the incision site is remote from and radial to the site of the expansion, the internal bladder of the expandable implant may be filled to tissue tolerance at the time of surgery.

If the incision site is not remote from, and radial to, the site of expansion, the wound should be stable before tissue expansion begins. However, a slight amount of inflation to fill the pocket space without tension to the tissue may be possible initially. Ensure that the expander/implant can be percutaneously filled and that there is no kink in the tube.

- a) Insert a new, sterile, 21-gauge needle (or smaller) into the injection site. Ideally, the needle should enter perpendicularly to the top of the injection site as shown in Figure 3.


FIGURE 3

- b) Penetrate the injection site until the needle is stopped by the needle guard, as shown in Figure 4.


FIGURE 4

- c) Air Aspiration: Initially remove as much air as possible. Fill the implant with 50cm³ sterile saline for injection, and then remove all remaining air.
- d) Fill the implant only with sterile saline for injection and only through the injection site. Fill carefully and only to patient and tissue tolerance. DO NOT underfill or overfill the implant beyond the range specified. Filling is typically performed at weekly intervals. A Patient Fill Volume Record card is provided with each implant for recording fill volumes and monitoring the expansion process (see back of this book).
NOTE: The suggested total fill volume is listed on the product labelling.

4. Removal of Remote Injection Site

After satisfactory tissue expansion, the remote injection port may be removed and the remaining fill tube sealed using the Fill Tube Plug Kit, available separately.

SALINE IMPLANTS

TECHNIQUE FOR USING IMPLANTS WITH DIAPHRAGM VALVE

1. Fill tube insertion

Holding the implant, push the secondary strap closure to one side of the valve entrance. Figure 5 shows the secondary closure in place and diaphragm valve seal closed.


FIGURE 5

Wet the tip of the fill tube in sterile saline for injection. Insert the fill tube by gently pushing the tip into the valve entrance as far as the tube flange permits. Figure 6 shows the diaphragm valve held open by the fill tube, for filling or emptying.

When the flange on the fill tube meets the periphery of the valve with a snap, the tube is in position for filling.


FIGURE 6

2. Air aspiration

After the fill tube is properly inserted, remove any air from the implant by aspiration with an empty sterile syringe attached to the fill tube.

3. Placement

Ensure incision is sufficiently large (particularly for textured implants) to facilitate insertion and avoid damage to the device. Inadequate pocket dissection increases the risk of deflation and implant malposition.

4. Filling

Use aseptic technique when adding/removing sterile saline, to prevent risk of contamination. Using a syringe filled with sterile saline for injection, inflate the implant. Fill only with sterile saline for injection, and fill to a volume within the recommended fill range specified on the product package labelling. DO NOT underfill or overfill the implant beyond the range specified.

5. Residual Air

After filling is completed, aspirate any residual air bubbles. Then use gentle traction to remove the fill tube from the valve, taking care to avoid damage to the shell or valve.

6. Diaphragm Valve Closure

Verify that the diaphragm valve is clear of particulates, and sealed. To help retard tissue ingrowth or fluid accumulation in the valve entrance, engage the strap closure as follows: using the thumb and forefinger, compress the valve seat and the strap to snap the valve plug into place as shown in Figure 5.

NOTE: The strap closure protects the valve after implantation, and the valve must be properly closed first to prevent leakage.

NATRELLE® 133 TISSUE EXPANDERS TECHNIQUE FOR USING NATRELLE® TISSUE EXPANDERS WITH MAGNA-SITE™ INTEGRATED INJECTION SITES

Natrelle® 133 tissue expanders are temporary devices, and are not intended to be used for permanent implantation. Tissue expanders should be removed once adequate tissue has developed. Tissue expansion in breast reconstruction typically requires four to six months. The total expansion period will vary depending on patient tolerance and desired flap size.

1. Natrelle® Tissue Expander Placement

Ensure incision is sufficiently large to facilitate insertion and avoid damage to the device. Inadequate pocket dissection increases the risk of deflation and tissue expander malposition. If the pocket is too small, the tissue expander may not have adequate room to unfold, increasing the risk of tissue erosion. Place the tissue expander flat end in the exact desired location, as the BIOCELL™ textured surface promotes tissue adherence to help immobilise the device. The MAGNA-SITE™ should be situated anteriorly, adjacent to the skin surface.

2. MAGNA-SITE™ Location

The MAGNA-FINDER™ is designed to assist with the location of the MAGNA-SITE™ injection site. While the injection site can be generally identified by palpation, always verify the location and orientation of the injection site with the MAGNA-FINDER™, as described below, before each filling.

- DO NOT store or use the MAGNA-FINDER™ near any loose metal particles, as they may attach themselves to the magnet.
- Be sure the magnet inside the MAGNA-FINDER™ moves freely without obstruction.
- Place the MAGNA-FINDER™ on the tissue overlying the implanted MAGNA-SITE™ contained in the Natrelle® 133 tissue expander.
- Following the direction of the arrow on the device, slowly move the base of the MAGNA-FINDER™ back and forth over the tissue surface, as shown in Figure 7. The magnet of the MAGNA-FINDER™ will indicate the location of the injection site on one plane.


FIGURE 7

- When the magnet on the MAGNA-FINDER™ is perpendicular to the base of the MAGNA-FINDER™, mark the tissue with a surgical marker at the two notches on the base of the MAGNA-FINDER™. The line between these two marks indicates location of one plane, as shown in Figure 8. The MAGNA-SITE™ will be at one point along this plane.


FIGURE 8

- To determine the exact location of the MAGNA-SITE™, turn the MAGNA-FINDER™ 90° and again, slowly move the MAGNA-FINDER™ back and forth. When the magnet locator is again perpendicular, make two more marks, as shown in Figure 9.


FIGURE 9

- Repeat this location process on several different axes to ensure the precise location of the injection site. The point at which all the marked lines intersect is the location of the MAGNA-SITE™. DO NOT re-sterilise or re-use the MAGNA-FINDER™.

3. Natrelle® Tissue Expander Filling

If the incision site is remote from the site of expansion, the tissue expander may be filled to tissue tolerance at the time of surgery. If the incision site is not remote from the site of expansion, the wound should be stable before tissue expansion begins. However, a slight amount of inflation to fill the pocket space without tension to the tissue may be possible initially.

- Insert a new, sterile 21-gauge (or smaller) standard 12° bevel hypodermic needle into injection site. Ideally, the needle should enter perpendicular to the top of the injection site as shown in Figure 10.

DO NOT penetrate the injection site at an angle.


FIGURE 10

- Penetrate the injection site until the needle is stopped by the needle guard.
- Fill the tissue expander only with sterile saline for injection, and only through the injection site after precise location of the MAGNA-SITE™ with the MAGNA-FINDER™. Fill carefully to a volume within the recommended fill range specified on the device package labelling and only to patient and tissue tolerance. The suggested fill volume is located on the patch of the tissue expander and is also listed on the device data sheet. Filling is typically performed at weekly intervals. A Natrelle® Expander Fill Volume Record is provided with each device for recording fill volumes and monitoring the expansion process (see back of this book).

INSTRUCTIONS AND PRECAUTIONS FOR REMOVAL

In case that removal of the product is necessary, the device should be removed from the patient without harming or breaking the outer shell of the product, if possible. If breakage of the outer shell occurs in an expander or implant filled with saline solution, no further action needs to be taken. If the outer shell of a silicone gel-filled implant is harmed or broken, all remnants of the silicone gel filling that might have exited the device must be removed from the patient's body mechanically or washed out. The wound and former implant pocket should then

be thoroughly rinsed with physiological saline solution. All devices removed from a patient's body should be treated as potentially contaminated, either safely disposed of or processed, and returned to Allergan as per the instructions below.

METHOD FOR REMOVING RUPTURED GEL FROM THE SURGICAL POCKET

In the event of implant rupture, the following technique is useful for removal of the gel mass. Wearing double talc-free surgical gloves on one hand, use the index finger to penetrate the gel mass. With the other hand, exert pressure on the breast to facilitate manipulation of the gel mass into the double-gloved hand. Once the gel is in hand, pull the outer glove over the gel mass and remove. To remove any residual gel, blot the surgical pocket with gauze sponges. Avoid contact between surgical instruments and the gel. If contact occurs, use isopropyl alcohol to remove the gel from the instruments. Ruptured implants must be reported and returned to Allergan as per the instructions below.

RETURNED GOODS POLICY

Product returns and exchanges must be authorised through your Allergan representative. Exchange value is based on time limitations. All package seals must be intact to be eligible for return or exchange. Returned products may be subject to a restocking charge. Certain products are non-returnable. For more information, please contact your Allergan representative.

REPORTING AND RETURN OF EXPLANTED DEVICES

Explanted devices associated with a complaint or injury must be reported and returned to your local Allergan representative with a Product Field Note (PFN). In order for the explanted product to be returned to the manufacturer it must first be decontaminated and a Decontamination Certificate completed and returned with the explanted device within an explant return kit. Do not return device if the patient has HIV or hepatitis or is known or suspected to carry another infectious agent. Please notify your Allergan representative in these circumstances.

PRE-DISINFECTING INSTRUCTIONS:

Gel-filled breast implants: Do not puncture the device.

For saline-filled, gel/saline double lumen breast implants and tissue expanders

- Remove any remaining saline solution from the device and vent the device by puncturing the shell with a needle, or nicking the shell with a scalpel. For gel-saline devices, puncture the saline-filled lumen, only.
- Mark the puncture by circling the site with methylene blue, skin marking pencil, or permanent marker.

Double lumen gel/saline implants that have been plugged with a Fill Tube Plug Kit

- Cut the tube approx. 1 cm from the Fill Tube Plug (implant side of the plug) and remove saline from the device, return both plugged tubing and the device.

DISINFECTING METHODS:

Autoclave:

- DO NOT use a prevacuum autoclave or ethylene oxide sterilizer.
- Set autoclave on the "slow exhaust" or "liquid" setting.
- Autoclave by the following gravity displacement cycle: minimum of 70 minutes at 121°C, 1kg/cm² (250°F, 15 psi).
- Open door slowly once the cycle is completed and allow device to cool to room temperature before preparing for shipment.

Bleach (use only if autoclave is not available):

Note: Do not use alcohol, *Cidex*® (glutaraldehyde), formaldehyde or other solutions for disinfecting.

- Mix one part of household bleach (10% sodium hypochlorite) with 9 parts of water.
- Completely submerge the explanted device in the solution for 60-120 minutes.
- Rinse thoroughly with water and dry the device.

A Product Field Note, shipping instructions, and a decontamination certificate within an Allergan return kit is required prior to dispatch of each explanted device. These can be obtained by contacting your local Allergan office/local distributor.

PRODUCT REPLACEMENT POLICY

In the event of non-iatrogenic loss of breast implant shell integrity within ten years from the date of implantation, the device will be replaced with an identical or equivalent device. To receive a replacement device, a Product Field Note (PFN) must be completed and submitted to your Allergan representative. Explanted devices must be returned in accordance with section above "Reporting and Return of explanted devices". The product replacement policy does not cover surgical or other expenses related to rupture, deflation, cosmetic revision, capsular contracture, or other adverse events. Contact your local Allergan representation for specific warranty details for your region.

LIMITED WARRANTY, LIMITATION OF LIABILITY, AND DISCLAIMER OF OTHER WARRANTIES

Allergan warrants that reasonable care was used in the manufacture and production of this product. Allergan has no control over the conditions of use, patient selection, surgical procedure, post-surgical stresses, or handling of the device after it leaves our possession. Allergan does not warrant either a good effect or against an ill effect following its use. Allergan shall not be responsible for any incidental or consequential loss, damage or expenses directly or indirectly arising from use of this product. Allergan's sole responsibility in the event that Allergan determines the product was defective when shipped by Allergan, shall be replacement of the product. This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether express or implied by operation of law, or otherwise, including, however not limited to, any implied warranties of merchantability or fitness for use.

INFORMED CONSENT

An Informed Consent Form is provided (see back of book). Please ensure that the patient receives the information from the section "INFORMATION THAT SHOULD BE PROVIDED TO THE PATIENT" and understands the information provided. The patient must realise that the surgical and post-surgical risks associated with implants cannot be completely predicted, even with the best medical manufacturing, technology and surgical care, and accept these conditions and limitations. Patients must fully inform their physician of their medical history, including any and all conditions that would contraindicate implant surgery or tissue expansion. Failure to inform their physician could result in significant surgical and post-surgical complications. The patient must decide herself on whether the expected benefits outweigh the said risks. If the patient decides that the expected benefits of the proposed implant surgery outweigh the risks, then she must take full responsibility for her choice to proceed with implant surgery. The two-part form located at the back of this book should be completed and signed once the patient has decided to proceed with implant surgery. This form allows both the patient and the surgeon to retain copies for their records.


ID CARD

As part of device tracking, Allergan is providing a patient ID Card. The information on the inside of this card (when folded) is specific to the device(s) the patient received (patients should keep this card for their records and carry it at all times to facilitate medical care in case of emergency). If a device is replaced, another card will be provided for the new device. The ID Card is located at the back of this book.

Instructions for the Surgeon:

Place one label from each product in the appropriate space on the ID Card (L or R). These labels can be found attached to the bottom of the main inner label on top of the inner product packaging. In the absence of a Patient Record Label, copy the Product code, REF, SN and LOT from the package label to the space provided. Fill in all the remaining sections. Give this entire document to the patient for her records. See "Graphical Symbols" for explanation of symbols.

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EN Natrelle® Informed Consent Form

- I confirm that I have read the *Natrelle®* implant/tissue expander patient information brochure.
- I confirm that I have had ample opportunity to discuss my concerns and options with my surgeon.
- I confirm that my surgeon has given me time to consider *Natrelle®* implant/tissue expander surgery.
- I fully understand the effects and risks associated with having *Natrelle®* implant/tissue expander surgery.
- I consent to undergoing *Natrelle®* implant/tissue expander surgery as discussed with my surgeon.

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